

The best choice for a perfect beauty care company is YOU & ME KOREA













NATURE CORRECTION MADE IN KOREA


Intended for facial tissue augmentation



Considering your beauty



A new generation of Cross-linked Hyaluronic Acid fillers DERMARCANE is aseptic products, highly tolerated, viscoelastic filler made with cross-linked hyaluronic acid, ideal for subcutaneous applications.





Best choice of total beauty care company is You and Me Korea Co., Ltd.

What is DERMARCANE dermal filler?

A new generation of Cross-linked Hyaluronic Acid fillers DERMARCANE is aseptic products, highly tolerated, viscoelastic filler made with cross-linked hyaluronic acid, ideal for subcutaneous applications.

Monophasic Cross-linked Hyaluronic Acid fillers DERMARCANE is safely injectable, safty, biodegradable medical device made with hyaluronic acid created from non-animal sourced bacterial biofermentation. The extremely low concentration of impurity and the bacterial endotoxin rate < 0.05 EU/ml minimise any hypersensitivity reaction.

It is a colourless monophasic gel contained in a disposable prefilled, ergonomic syringe with a Luer Lock adapter that makes it easy to handle.

Minimise any hypersensitivity reaction DERMARCANE is manufactured under the strict controls of the certificates ISO 13485 and Ministry of Food and Drug Satety by Korea.

DERMARCANE has an immediate lifting effect, thanks to it's intrinsic viscoelastic and hydrophilic properties. The treatment result is very natural and immediately visible. 24mg/ml concentration Hyaluronic Acid allows reducing any type of deep wrinkles, reshaping the face contour and lips, and hydrating the skin in depth.



Formulation



Filling



Stoppering



Sterilization



Visual inspection



Labelling



Blister packaging



Quality control



DERMARCANE Products

Intended for facial tissue augmentation







DERMARCANE Grey with lidocaine

Volume creation of the facial contours, forehead, chin, cheeks, nose bridge and etc.



DERMARCANE Red with lidocaine

Tissue augmentation of the deep facial lines and skin depressions, lip, nose, nasolabial folds, marionette lines and etc.



DERMARCANE White with lidocaine

Line remediation of the medium facial lines and skin depression, glabellar lines, lip contour and etc.



DERMARCANE Green with lidocaine

Wrinkle correction of the superficial facial lines and skin depressions, perioral lines, forehead finelines and etc.

	Grey	Red	White	Green
Indications	Volume creation and facial contours	Deep facial wrinkles and folds	Medium facial wrinkles	Superficial facial wrinkles
Depth of injection	Subcutaneous	Deep dermis	Mid-dermis	Superficial dermis
HA concentration	24mg/ml	24mg/ml	24mg/ml	24mg/ml
Lidocaine concentration	3mg/ml	3mg/ml	3mg/ml	3mg/ml
Fill volume	1.1ml	1.1ml	1.1ml	1.1ml
Needles	2 x 25G ^{1/2'}	2 x 26G ^{1/2'}	2 x 27G ^{1/2'}	2 x 30G ^{1/2'}
Self life	24months	24months	24months	24months



DERMARCANE Feature



* EDQM (European Directorate for the Quality of Medicines and Healthcare)

Frequently Asked Questions

- [Q] How often will I need a treatment?
- [A] People tend to have top-up treatments once they notice the benefits beginning to wear off, usually anytime between 6 and 24 months.
- [Q] Who is suitable for Hyaluronic Acid Dermal Fillers?
- [A] Almost everyone, providing they're not pregnant or breastfeeding.
- [Q] Are there any long-term problems?
- [A] No. One of the real beauties of Hyaluronic Acid Dermal Fillers is there aren't any long-term problems because it's a natural treatment that's safely broken down by the body over time.
- [Q] Are there any side effects?
- [A] You might experience some minor swelling or redness which should last no more than a day or two, or very occasionally up to one week. A tiny proportion of people experience some bruising around the injection site.
- [Q] Are the injections painful?
- [A] The injections aren't painful, because our Hyaluronic Acid Dermal Fillers include anaesthetic called "Lidocaine".